

# UV-Visible Spectroscopy for concentration calculation

Date: 2020-01-21

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

1 / 4

Goal : Establish a calibration curve of  $\text{Ti}_3\text{C}_2$  MXene for UV-Visible spectroscopy

Procedure :

- Perform serial dilution of bulk  $\text{Ti}_3\text{C}_2$  MXene suspension from latest synthesis (see [\[Experiment\] MXene synthesis IV](#)) at  $13.8 \pm 0.3$  wt% with a target minimum volume of 5 mL per sample
  - Use disposable 1 mL plastic pipette to transfer a volume of one suspension to an empty glass vial
  - Record mass of suspension added
  - Add a quantity of deionised water to achieve desirable concentration with a separate pipette and record total mass
- Ensure homogeneity with manual shaking
- Add homogeneous solution to quartz cuvette and insert into spectrophotometer ensuring correct orientation and dry outer faces
- Perform scan with speed 480 nm/min and 1 nm data intervals in 200-700 nm range

Results :

Sample N° / #	Mass of previous suspension added / g	Mass of deionised water added / g	Concentration / wt%
1	0.0600	7.2891	$1.127 \times 10^{-1}$
2	0.3397	6.6106	$5.507 \times 10^{-3}$
3	1.1368	3.7625	$1.278 \times 10^{-3}$
4	2.4383	4.0706	$4.786 \times 10^{-4}$
5	3.7350	3.8273	$2.364 \times 10^{-4}$
6	4.3455	3.1185	$1.376 \times 10^{-4}$
7	0.4102	6.4425	$8.239 \times 10^{-6}$

.csv files are xy data of wavelength vs absorbance, .sp and .rtf files are generated

# UV-Visible Spectroscopy for concentration calculation

Date: 2020-01-21

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

2 / 4

by Perkin-Elmer UV Win Lab software and are spectra and a report file, respectively. .png files are plots of the data.

## Conclusions:

Samples of concentration  $\approx 5.507 \times 10^{-3}$  wt% are not suitable for this spectrophotometer, as the absorbance exceeds the limits of detection. In samples of concentration  $\approx 1.278 \times 10^{-3}$  wt%, peaks are present across all datasets at 322 and 562 nm. Calibration curves of concentration versus absorbance at each of these wavelengths gives datasets which are well-represented by straight-line plots when excluding the unsuitable samples.

## Attached files

21-01-20-full-report-rtf.rtf

sha256: 389c76cc80df924eea5341b6d44c1489bb2031618a30caebfbac651dc3598245

MXene-conc.-3.Sample.Raw.csv

sha256: 3a9b62006703e6ce00685b3d651254c158a2b6ec83475ea1685a0f021156c926

MXene-conc.-3.Sample.Raw.sp

sha256: 417ca301b3928474447afcd03d7f2d2d496af721a44b6e2bfa1869acc6d83f01

MXene-conc.-4.Sample.Raw.csv

sha256: b0b422ffb415c76fe070a0f906748c0c559ac944250f4a1383c2d24d351eb78c

MXene-conc.-4.Sample.Raw.sp

sha256: 263ca0a62150412bb3d37cc0714e983872d5349b8a7a243afcee280694b770d8

Mxene-conc.-5.Sample.Raw.csv

sha256: a03ec8c3aa9654594b439ef569df60a5e5a9740163e3b7c39d51bac519e97946

Mxene-conc.-5.Sample.Raw.sp

sha256: f07f91cc80b79e19712464b35f8db60ba41a7f8852751580ced1db0b50618aba

# UV-Visible Spectroscopy for concentration calculation

Date: 2020-01-21

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

3 / 4

Mxene-conc.-6.Sample.Raw.csv

sha256: 186e982554fca1b70d8f91940465cc6cd1bff8148046d39ccb33ba35cb2c7b2e

Mxene-conc.-6.Sample.Raw.sp

sha256: 044e2afc0c0a55ec7ecba81ad46bcb120ac4f19a15c6b42b01489f7ae669577f

MXene-conc.-7.Sample.Raw.csv

sha256: 87d54404286e180a1893c911acbf7ba9bbdd51d697681c8e5165ed97443d61d4

100-or-0-Absorbance-Baseline.Correction.Raw.csv

sha256: 8154e252dac84523a939b844438a06580380da80e9a422d7ecd8c0e7af10e936

MXene-conc.-7.Sample.Raw.sp

sha256: bc04d441e86b313d1376a50cbcd4292713e154a8e398c50296b1fc2747fc22db

100-or-0-Absorbance-Baseline.Correction.Raw.sp

sha256: 1692402d1d3b518251c201e97b71364bee88a03dbbb9d23b1768f92ac7bf3180

MXene-conc.-1.Sample.Raw.csv

sha256: 814a4f1c8702b039a8fa7d7682cf23747cf082bc2c6ba3b2bb13cbc8646b2d1b

MXene-conc.-1.Sample.Raw.sp

sha256: 91f2ce78633a286450b72dd2eaaaa0930709c52c9deabc1f29c8548a0769865d

MXene-conc.-2.Sample.Raw.csv


sha256: cd5df7c9bd6dc64dcf4bf0b0e558b910e29eec4577ecccc0a71c384ffadf897c

MXene-conc.-2.Sample.Raw.sp

sha256: edc71c86dbce4267c4cb0ab824ea6e1455bc4ee802b64cdc778f787743430fb2

UVVis200121.png

sha256: c0b784a94173c53cbf5f3ab7e0ce8e76278288314f4eb927eda94b8d0bc40d15



UVVis200121ylim.png

sha256: f590b5b323edbd4465933f35825cfe29d7b73b9b96be7f26fedb687129c47dbf


# UV-Visible Spectroscopy for concentration calculation

Date: 2020-01-21

Tags: 14/01/2020Synth UV-Vis Concentration Calibration


Created by: James Bird

4 / 4




UVVis200121\_322nmAbs.png

sha256: 3bbe5d088b3afbf58d6bdb9b4047359dca91fd378909c736a3a399556b466a1d



UVVis200121\_562nmAbs.png

sha256: 5bb3c8bcec005b8ef5a91d147cf45a9304508e972c9fbf6049d8bd31e452921a



Unique eLabID: 20221014-e58da5574b40d7138708c359781ccf281e209b2b

Link: <https://frankel-elab.manchester.ac.uk/experiments.php?mode=view&id=37>