

# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05

Tags: *PSD DLS Article 1 Synthesis Optimisation*

Created by: James Bird

1 / 7

Goal : Obtain particle size distributions (PSDs) of  $\text{Ti}_3\text{C}_2$  MXene nanoparticles in aqueous suspension using DLS

Procedure :

## Sample preparation

- Bulk concentration  $\text{Ti}_3\text{C}_2$  product suspensions, of concentrations given in [\[Experiment\] Optimisation of MXene synthesis - Yield calculations](#), are diluted from to give a volume of ~1.2mL at 0.1 wt%. Dilutions are detailed in the table below, where suspension density is approximated to  $1 \text{ g cm}^{-3}$

Run N° / #	Volume of bulk suspension / $\mu\text{L}$	Volume of deionised water / $\mu\text{L}$
1	51	1149
2	98	1102
3	50	1150
4	1091	109
5	62	1138
6	51	1149
7	308	892
8	250	950
9	52	1148
10	308	892
11	500	700
12	156	1050
14	107	1093
15	194	1000
16	65	1135

- Dilution is carried out by extracting a quantity of the bulk suspension using a

# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05

Tags: PSD DLS Article 1 Synthesis Optimisation

Created by: James Bird

2 / 7

Gilson Pipetman P200 micropipette and D200 tips and transferring to a glass vial. The desired quantity of deionised water is measured in the same way.

- If a stated target volume exceeds 1 mL, then a 1mL disposable pipette is used to transfer that quantity
- Suspensions are homogenised throughout with vortex mixing
- Final concentration suspensions transferred to cuvette with Pasteur pipette - cuvette exterior gently dried if necessary

## DLS operation

Standard Operating Procedure (SOP) settings (size measurement type):

- Narrow band filter fitted
- Material properties (RI = 1.7, absorption = 0.9)
- Water dispersant ( $\eta$  = 0.8872 cP, RI = 1.330)
- Use dispersant viscosity as sample viscosity
- Temperature = 25 °C with 120 s equilibration time
- DTS1070 folded capillary cell *only*
- 173 ° backscatter measurement angle
- Automatic measurement duration
- Three measurements per sample
- Automatic attenuation selection and positioning method seeking optimum
- General purpose (normal resolution) analysis model *only*

Results :

All measurement outputs and variables indicative of data quality are detailed in the table below. All individual measurements returned a polydispersity index (PDI) > 0.1, which defines the threshold below which outputs can only be compared quantitatively. Despite this, the combination of three measurements for one

# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05

Tags: PSD DLS Article 1 Synthesis Optimisation

Created by: James Bird

3 / 7

sample did meet the data quality criteria defined in the analysis program (Zetasizer) : all data quality indicators are highlighted in green in that row (Run #7).

Run N° / #	Z-average / d.nm	Z-average std / d.nm	Number mean / d.nm	Number mean std	PDI / dimensionless	PDI std	Derived mean count rate / kcps	Derived mean count rate std	Intercept	Intercept std	In range / %	In range std
1	3641	515.3	2051	425.1	0.480	0.266	9771.5	3579.9	0.952	0.021	77.1	1.64
2	2888	633.1	344.4	97.41	0.986	0.025	8416.3	1262.7	1.42	0.194	81.9	1.18
3	1779	338.6	938.4	123.9	0.406	0.211	4813.3	702.7	1.01	0.091	81.2	1.01
4	286.2	65.6	92.83	35.67	0.440	0.053	599.3	25.2	0.837	0.008	89.2	1.83
5	1750	314.3	852.5	212.9	0.599	0.088	2910.5	353.7	0.920	0.045	77.1	2.35
6	2208	268.7	463.7	147.7	0.945	0.017	5305.3	913.4	1.17	0.095	76	5.14
7	540.5	6.086	277.3	15.43	0.409	0.031	373.9	0.8	0.841	0.003	96.7	0.84
8	812.9	134.4	251.8	26.57	0.691	0.040	15138.3	1336.4	1.07	0.035	88.5	1.80
9	2952	199	773	227.1	0.377	0.011	277.5	8.9	0.769	0.01	91.3	2.5
10	902.9	195.3	198.1	34.78	0.82	0.162	109.6	20.7	0.968	0.062	72.6	2.86
11	4935	942.2	2419	233.1	0.574	0.378	4473.5	821.3	0.965	0.086	69.7	1.2
12	1728	251.8	685.6	173.5	0.737	0.169	6584.3	871.2	1.04	0.122	77.9	0.85
14	3264	261	1148	642.7	0.534	0.438	2058.7	1107.1	0.956	0.057	65.5	2.62
15	1628	248.3	626.7	66.04	0.789	0.133	1459	225.8	0.991	0.061	76.1	2.46
16	2591	201.2	859.2	328.7	0.743	0.235	640.4	293.2	0.913	0.053	67.4	4.34

Runs #4 & 9 are highlighted in yellow, as three of the four quality indicator variables are in a suitable range, and the fourth (percentage in range) is greater than all other run averages. Common to all of Runs 4, 7 and 9 in the intensity PSDs is a sharp drop-off in density at the larger end of the particle dimension scale: this corresponds to a poor cumulants and distribution fit to the correlogram on longer timescales, corresponding to those larger particle dimensions. Run #5 is also highlighted in yellow as three of the four quality indicators are suitable, and shares a similar issue, in that the cumulants analysis fits well, but distribution analysis does not on longer timescales.

.dts is the raw datafile which can be read into the Zetasizer software program, .csv is an exported, comma-separated summary of the .dts datafile and .png files are the plotted data; particle diameter (on a log scale) is plotted against the mean percentage of that diameter interval contributing to the intensity-based and

# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05

Tags: PSD DLS Article 1 Synthesis Optimisation

Created by: James Bird

4 / 7

number-based particle size distribution for each synthesis run. Vertical lines spanning the whole plot height are mean Z-average and number mean values (quoted above), where the regions of matching colour spanning left and right of this value correspond to its standard deviation.

## Conclusions:

The intercept regularly exceeding 1, PDI exceeding 0.6 and low in-range percentages suggest any one of a number of issues with the samples, including large or sedimenting particle presence, unsuitable suspension concentration or sample polydispersity too high for DLS measurements. In combination with the zeta-potential measurements (detailed in [\[Experiment\] Zeta-potential measurement of MXene suspensions II](#)), the first option can be ruled out for 11 of these measurements. Certainly, concentration optimisation should have been carried out for each individual synthesis run to ensure no impact stemmed from this variable, as this has a particle-size dependence on the optimal value. However, the fact that nearly all samples give a high multimodal fit error does suggest that this sample type may be inherently too polydisperse. The software does have an option to extend duration for large particles, which also should have been tested.

## Attached files

2021-05-06\_Size.csv


sha256: 428fc7bc5fab78c71aa9bb683d67a1a3d1ba36336afdf7b27f575eb03f8871c7

2021-05-06.dts

sha256: c4df992aebdc305b0d5c60de5a2e071480e3262416b0fc7de0ad346e1d448dbd

DLS\_SynthOpt\_PSDComp\_Run11.png

sha256: 83cc82919d77522f5dd4f723c88922cc117c024d6aa24a1a2ca429391f96da2f



# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05


Tags: *PSD DLS Article 1 Synthesis Optimisation*

Created by: James Bird

5 / 7


DLS\_SynthOpt\_PSDComp\_Run10.png

sha256: 9ccd4cba3f6faf99996cd49f42753e01e9f218bb57263ec08374d658fd7533cd




DLS\_SynthOpt\_PSDComp\_Run4.png

sha256: aaf9204de2781f671dd64430e1d63de0352551e3e9f491d7eeb6306d6f0d5c5d




DLS\_SynthOpt\_PSDComp\_Run9.png

sha256: 2d3f3d318ff967ca8991de38544b8baf1d35113ef5e7cb2b01b550036e2c2168




DLS\_SynthOpt\_PSDComp\_Run8.png

sha256: e356b8a80d0a9c11233e3e84dcca50652a1592ecfbb8134ba4fb9012e1bd56bc




DLS\_SynthOpt\_PSDComp\_Run7.png

sha256: 1965842b721a492375f5efc3f640b7fdc2ec27e57ca9ecf00375ed84373a19fd



DLS\_SynthOpt\_PSDComp\_Run6.png

sha256: 79c75bd0d666704073f26fe3837cd3614d35d21ca6ba7c7e797d58a311a3383e



# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05


Tags: PSD DLS Article 1 Synthesis Optimisation

Created by: James Bird

6 / 7


DLS\_SynthOpt\_PSDComp\_Run5.png

sha256: ae4f7c0de7ebbdd6db56c5db4499c6dd2ab9fa5ca135b19a594446ce9ed27330




DLS\_SynthOpt\_PSDComp\_Run3.png

sha256: 01f99004eea2c093b55c87aae72f8deaa31492d0907977261c5615132b1c3c8b




DLS\_SynthOpt\_PSDComp\_Run1.png

sha256: 26e5d7711a42123c824bef8a54450b93cf487d7b20345f4ff19e1ab407f7b70a




DLS\_SynthOpt\_PSDComp\_Run2.png

sha256: 6e67a2d607f418e5e8cd6344219ba519d713df4d34c2eceabff810dbd70c41f7




DLS\_SynthOpt\_PSDComp\_Run16.png

sha256: 91c9f94f167c93fe88d5ed573d7d551d7cc71cabb93a26be2db1dd51f48931bf



DLS\_SynthOpt\_PSDComp\_Run15.png

sha256: d2fa84c6f4cf5ae2b451a0dc4083f847fee579fc4b2ee63e3f388ae285922690



# Dynamic Light Scattering (DLS) for nanoparticle size-distribution acquisition III

Date: 2021-05-05


Tags: *PSD DLS Article 1 Synthesis Optimisation*

Created by: James Bird

7 / 7


DLS\_SynthOpt\_PSDComp\_Run14.png

sha256: 7776599bc6b5188a45457d5e5480b85c94305ad8e28aea805483e1c947557a1a



DLS\_SynthOpt\_PSDComp\_Run12.png

sha256: 6c26dc7e9c0ab4043c0c4a1a32b5938cdcb5859eaaaa32a101ddcc35194dee8d



Unique eLabID: 20230124-456a4f1c09da999803a6f0e7e5eef9c5ebcdc17

Link: <https://frankel-elab.manchester.ac.uk/experiments.php?mode=view&id=80>