



HARD X-RAY PHOTOELECTRON SPECTROSCOPY



The world's first high-throughput laboratory HAXPES at 9.25 keV

- Excillum gallium metal jet (Ga K α) X-ray source
- Bespoke monochromator
- Micro-focussed to 50 μm to gain 3 orders of magnitude greater flux compared to laboratory Al K α sources (see Regoutz *et al.*, Rev. Sci. Instr. **89** 073105 (2018))
- Combined with EW4000 high kinetic energy electron analyser
- Connected under ultra-high-vacuum to a high-throughput XPS with UV lamp, Al/Mg dual anode, monoatomic and cluster etching
- Heating and surface preparation chamber
- Vacuum suitcase for sample transport without exposure to atmosphere

Enabling chemical state information for thin films, buried interfaces and active layers, multi-layered structures, and bulk material information



- HAXPES applies the powerful tool of XPS up to depths of at least 50 nm below the surface
- Analysis of the inelastic background can enable depth-profiling through up to an estimated 20 x IMFP for depths over *hundreds of nm*
- Much less sensitive to surface roughness and contamination

A new tool for advanced materials researchers and industry

- Bringing the technique of HAXPES outside of national-level synchrotron radiation facilities opens up great possibilities for application to a wide range of research across multiple disciplines
- Expect a revolution in uptake as previously occurred with lab-based near-ambient pressure (NAP) XPS
- Establishment of the technique includes a 3-year collaboration between Scienta Omicron, The University of Manchester and the National Physical Laboratory (see data in the left-hand figure)
- Standardisation of HAXPES – establishing sensitivity factors, sampling depths, and analysis protocols for a wide range of materials systems

Access also available through the XPS National Research Facility, www.harwellxps.uk






1s core level				
Element	BE (eV)	HAXPES	XPS	Sampling depth (nm) est.
Li	55	66.9	13.5	(3 x IMFP)
Be	112	45.9	9.3	
B	188	45.6	9	
C	285	51.9	9.9	$h\nu$
N	410	51.6	8.7	HAXPES 9252 eV
O	540	50.1	7.8	XPS 1486 eV
F	695	48.9	6.9	
Ne	870	49.2	5.7	
2p core level				
Na	1070	74.1	6.3 BE (eV)	
Mg	1303	47.7	2.4	
Al	1559	37.8		49 54.3 11.4
Si	1839	37.2		72 44.1 9.3
P	2145	37.5		99 44.4 9.3
S	2472	34.8		136 46.5 9.3
Cl	2822	35.1		164 44.4 9
Ar	3205	34.2		200 47.1 9.3
K	3608	69.6		248 48 9
Ca	4038	40.5		294 103.8 18.9
Sc	4492	27.3		346 64.2 11.4
Ti	4966	20.7		398 46.2 8.1
V	5465	16.2		453 37.8 6.6
Cr	5989	13.2		512 32.4 5.4
Mn	6539	11.1		574 29.4 4.8
Fe	7112	8.7		638 28.5 4.5
Co	7709	6.3		706 27.3 4.2
Ni	8333	4.5		778 25.5 3.6
Cu	8979	1.8		852 25.5 3.3
				932 25.2 3

HAXPES sampling depths

Using photoelectron peaks (3 x IMFP) compared to Al K α XPS. HAXPES enables measurement of deeper core levels, and for larger atoms, using different core levels gains information from different depths.

Standardisation of HAXPES: sensitivity factors, sampling depths, buried layers

- Survey spectrum ($h\nu = 9$ keV) for $\text{Ir(ppy)}_2(\text{acac})$ containing 3% iridium; deeper core levels (3 spd) have very large sensitivity factors. Auger transition peaks are at much lower kinetic energies w.r.t. the photoelectron peaks. *Theoretical RSFs obtain the correct atomic percentages.*
- Ir 4spdf with O, N, C 1s photoelectron peaks.
- Ir 3spd region for a silicon wafer with 200 nm organic underlayer (Irganox 1010), then with an 18nm layer of Ir complex, then with various organic overlayers up to 140 nm (as shown in the drawing). Core levels are lost for overlayers > 50 nm thick, however, *the inelastic background associated with Ir 3spd electrons still contains information which can be modelled.*
- h-h Ir 3d region for various organic overlayers (140 nm, 100 nm, 75 nm, 50 nm, 25 nm as labelled). The overlayer thickness is obtained by modelling (yellow lines) the inelastic background associated with Ir 3d, using Sven Tougaard analysis (Surf. and Inter. Anal. **11** 453 (1988); Journal of Surface Analysis **24** 107 (2017)) – QUASES software (www.quases.com).



Manufacturer: Scienta Omicron GmbH

Academic Lead: Prof Wendy Flavell

Technical Lead: Dr Ben Spencer

Research & Facilities Manager: Dr Helen Ryder

Location: G318 Photon Science Institute, Alan Turing Building, The University of Manchester

“Due to its relatively unlimited electron escape depths, HAXPES has emerged as a powerful tool that has general application to the study of the true bulk and buried interface properties of complex materials systems. Its areas of application are thus growing exponentially compared to more traditional measurements at lower photon energies.”

Joseph C. Woicik, NIST

