

UV-Visible Spectroscopy for concentration calculation I

Date: 2022-01-31

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

1 / 5

Goal : Re-establish a calibration curve of Ti_3C_2 MXene for UV-Visible spectroscopy by performing scans across the complete wavelength range of the spectrophotometer

Procedure :

- Produce an intermediate concentration suspension between the target values and the bulk Ti_3C_2 MXene suspension from latest synthesis (see [\[Experiment MXene synthesis IV\]](#)) up to ~ 25 mL
 - Use disposable 1 mL plastic pipette to transfer a volume of one suspension to an empty glass vial
 - Record mass of suspension added
 - Add a quantity of deionised water to achieve desirable concentration with a separate pipette and record total mass
- Use above method to produce all more dilute suspensions from this one intermediate (Sample #13 below) with a target volume of ~ 15mL per concentration
- Ensure homogeneity with manual shaking
- Add homogeneous solution to quartz cuvette and insert into spectrophotometer ensuring correct orientation and dry outer faces
- Perform scan with speed 480 nm/min and 1 nm data intervals in 200-1000 nm range
- When samples not in use they are degassed, sealed and stored refrigerated

Results :

Sample N° / #	Mass of suspension added / mg	Mass of deionised water added / g	Concentration / wt%
13	102.1 (from bulk suspension)	25.9496	5.430×10^{-2}
8	183.3	16.3722	6.079×10^{-4}
9	222.3	16.3410	7.387×10^{-4}

UV-Visible Spectroscopy for concentration calculation I

Date: 2022-01-31

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

2 / 5

10	264.0	16.7162	8.576×10^{-4}
11	286.4	15.4176	1.009×10^{-3}
12	337.4	16.0314	1.047×10^{-3}
3.1	343.1	14.6003	1.276×10^{-3}
4.1	125.4	14.1804	4.802×10^{-4}
5.1	80.7	18.2449	2.402×10^{-4}
6.1	29.7	11.5107	1.401×10^{-4}
7.1	35.5	31.7034	6.080×10^{-5}

.csv files are xy data of wavelength vs absorbance, .sp and .rtf files are generated by Perkin-Elmer UV Win Lab software and are spectra and a report file, respectively. .png files are plots of the data.

Conclusions:

All sample concentrations are suitable for this spectrophotometer. Consistent peaks are again identified present at 322 and 562 nm. On extending the data range to 1000 nm, another consistent peak is identified at average position 742 nm.

Attached files

KaiKai-Ti3-MILD.rtf

sha256: 64344ed4c7e28f00953e465d8fa5d03dd2bcfb44a49853af9a1911ac419b1068

8.Sample.Raw.csv

sha256: fcc99f7a865babe0fb98eeb94eb8f63d0840307bef291c5e6928db03fe9360a8

9.Sample.Raw.csv

sha256: 9182dc4b28f7bb68119a492444c42854f8aa949cfee1ee975356869916869ff3

UV-Visible Spectroscopy for concentration calculation I

Date: 2022-01-31

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

3 / 5

10.Sample.Raw.csv

sha256: 3120a0da672ab4d3abec92bf323e477e10ffbb9fe77fac5ca6aa85a4a7e0cee5

11.Sample.Raw.csv

sha256: 5d2d511fde649f54af7410e28f599dce2828e0644833e0baf3e8cdb648881371

12.Sample.Raw.csv

sha256: 8f2e22ab6c58e8f4dfcabd5f2f55c3f4d5e046c0a88989662415089378b6506f

100-or-0-Absorbance-Baseline.Correction.Raw.sp

sha256: cbd12be55fc80dd9683159bad457ab4d9d64883a9a678438c3c76434390ba1cc

7.1.Sample.Raw.sp

sha256: 1df0c2d29a2695867accd8d34c010881dc54919dec5dfb083375253a29457382

3.1.Sample.Raw.sp

sha256: 11a7cd1b77edcf859d0545c369958c8f3089644468c5b6d20ce8f51436b446f0

4.1.Sample.Raw.sp

sha256: 4f21ebd38a4ff38f89da88c8008a1e2a58c85afde10ba11af54d0094759d9101

5.1.Sample.Raw.sp

sha256: ee006a3ef950b131164dff85505ef5e1b2acd2c373b0c66c1eec73ea2134e30

6.1.Sample.Raw.sp

sha256: b7f61a91d6b528101a3bdc3db78a7b64358228816adf6b8166bfc32dfa0cf8e8

10.Sample.Raw.sp

sha256: 42fabbb72ca4b72e2c6b37ccf5d74e9d59b045d5d40c22394f18672d65f8c9ef0

11.Sample.Raw.sp

sha256: e26611f231902b5d88ba9d3e65c580f84b1aa221a1f65a53b99f8ecbe5aff06d

12.Sample.Raw.sp

sha256: a7c967e909adc1a3f7d80795154047e0c8fc0a6f728670e3a5b96c56c2bbe3b0

8.Sample.Raw.sp

sha256: 65afd26179d957983b1bf80950bac9f3f13f26b3e5a9495472d2f3d7f8ee2e9c

UV-Visible Spectroscopy for concentration calculation I

Date: 2022-01-31

Tags: 14/01/2020Synth UV-Vis Concentration Calibration

Created by: James Bird

4 / 5

9.Sample.Raw.sp

sha256: d42b4acec0e3e56865c40d38ae9cf75a23bfe79c21a54d6a25d44bc7890f7f45

100-or-0-Absorbance-Baseline.Correction.Raw.csv

sha256: 0178d80b4aec7bbb5eef0a17a9834b201bf899b41499200fdfff8b68c7e6d866

3.1.Sample.Raw.csv

sha256: 07811b46a6ed6cf4430985cdfa3d86ecc0624803948879f3ea9e6c5581d114e4

4.1.Sample.Raw.csv

sha256: 88e5fc33c435503c19beba782a25541e15ec75f820fb929258434601aa556ef7

5.1.Sample.Raw.csv

sha256: 3e362e2d04c0e1741f46d1c57fb57bf8a3883ff797e9e4b5e8af5824135c5131

6.1.Sample.Raw.csv


sha256: 999adb13436610cf94bb42ff9d3f689e84a0d527494347919f7c2dcdafbe2fa0

7.1.Sample.Raw.csv

sha256: 4b746d0338e0aecf7b3e93588923179373b97021f39bbe3bbbe43970aad932cc


UVVis200204_562nmAbs.png

sha256: 14159aa466c4eb0f443c91d080bfa50d573cf66ba3110eb59cde977952eb459c



UVVis200204_742nmAbs.png

sha256: 5e7a7bb536aa356e3bbd6ca637e48371922fa8c60c056d790ae1424df133aa9d



UVVis200204.png

sha256: 134896f565b65b62e005d3f70fe2b05fb4c3a017789a672f0782bb438cc939a4


UV-Visible Spectroscopy for concentration calculation I

Date: 2022-01-31

Tags: 14/01/2020Synth UV-Vis Concentration Calibration


Created by: James Bird

5 / 5



UVVis200204_322nmAbs.png

sha256: 15d02113dea84475de949d7157dd1b4e469d0f121ed70717b09bbacff653acf1



Unique eLabID: 20221017-5c8821f5113c2fec062fd6feaca89df2c7403eb5

Link: <https://frankel-elab.manchester.ac.uk/experiments.php?mode=view&id=43>